


# Tonoclear


STERILE



Part Number: 2414-P-7045



*Tonoclear mires image*



*Competitor mires*

## Product Overview:

A single-use tonometer prism to be used in combination with a GAT, to determine the IOP. The Tonoclear prism is designed to accommodate variations in light intensity from the slit lamp to ensure a clear view of the mires at all times. The individually packed prisms come pre-assembled to reduce any risk of a cross-contamination and hassle in your workflow.

## Features:

- Single-use prisms (box of 100)
- Pre-assembled
- Clear visualisation of the mires
- Attenuating sleeve to reduce sensitivity to light variations
- Individually packaged and sterilised with Ethylene Oxide
- Sterilisation validity for 3 years
- Compatible with all GATs

For more information, visit [www.keeler.co.uk](http://www.keeler.co.uk) or email us at [info@keeler.co.uk](mailto:info@keeler.co.uk)